



# RNA EXTRACTION FROM LUNGS

**Tanja Thurn: Apr 25<sup>th</sup> 2011**  
Bowdish Lab, McMaster University  
Hamilton, ON, Canada  
[www.bowdish.ca](http://www.bowdish.ca)

## BACKGROUND

- For isolation of RNA from tissues with Trizol reagent cat# 15596-018, Invitrogen

## NOTES

- Perform this procedure in a fume hood, wear goggles, gloves and coat.

## EQUIPMENT

- Equipment:
  - o Fume Hood in MDCL 4077
  - o Bullet Blender for homogenization
  - o Vortex
  - o Centrifuge
  - o Heat-block
  - o Tubes for Bullet Blender and Beads (3.2 or 1. mm beads)
  - o Pipettes and Tips, Tubes
  - o illustra RNAspin Mini kit from GE Healthcare – optional for cleanup and DNase treatment
- Reagents:
  - o RNA later
  - o Trizol reagent
  - o Chloroform
  - o Isopropanol
  - o 75 % Ethanol (prepared with RNase free water)
  - o RNase free water

## PROTOCOL

(adapted from Invitrogen [http://tools.invitrogen.com/content/sfs/manuals/trizol\\_reagent.pdf](http://tools.invitrogen.com/content/sfs/manuals/trizol_reagent.pdf))

- o Harvest lungs and put in **RNA later overnight at 4°C**, can be frozen after that at -80°C or flash freeze on dry ice – cells must be lysed in order to proceed with Trizol
- o Thaw, put lung in **1 ml of Trizol** (per 50-100 mg of tissue)
- o Add 5-8 beads and close tube very tightly!!
- o **Homogenize** in Bullet Blender for 5 min (option to freeze at this point)
- o Take suspension and incubate homogenized sample for around 2 more min at room temp. to permit complete dissociation of nucleoproteins

- **Add 200 ul of chloroform** per 1 ml of Trizol reagent (if frozen proceed with Trizol extraction protocol as soon as thawed by adding chloroform)
  - Vortex and incubate at room temp. for 2-3 min.
  - Centrifuge at no more than 12.000 x g for 15 min at 2-8°C
  - Take aqueous phase that contains RNA and put in new 1.5 ml tube
  - Precipitate RNA by adding **500 ul of Isopropyl alcohol** per 1 ml of Trizol reagent
  - Incubate sample for **1 hour at -20°C** or for **10 min. at room temp.**
  - Centrifuge at no more than 12.000 x g for 10 min at 2-8°C
  - Remove supernatant and wash RNA pellet with **1 ml of 75 % ethanol** (prepared with **RNase free water**)
  - Vortex and centrifuge at no more than 7.500 x g for 5 min at 2-8°C
  - Briefly air-dry pellet, do not over-dry and dissolve RNA in RNase free water, incubate at 55-60°C for 10 min
- I usually resuspend it in 90 ul

### Optional:

#### Clean-up and DNase treatment on column illustra RNAspin Mini kit from GE Healthcare

- Add 3.5 volumes of **buffer RA1** per 1 volume of sample, vortex (315 ul if pellet resuspended in 90 ul)
- Add **ethanol (95-100 %)** at 3.5 times the volume of sample, vortex
- Load sample directly on the **blue column** (RNAspin Mini column, GE Healthcare) and follow protocol as described in the booklet from the kit including digestion with DNase and at the elution step adding RNase inhibitor.
- Elute the RNA in 60-70 ul of pre-warmed (65 °C) RNase free water. First elution with 30 ul and second elution with 30-40 ul
- Spec RNA
- **Prepare your cDNA:** if not done on same day - **aliquot sample** (avoid freeze-thaw cycles with RNA) for convenient amount for cDNA prep: maximum 2 ug total RNA for cDNA prep in preferably 10 ul total volume

### LINKS AND REFERENCES

[http://tools.invitrogen.com/content/sfs/manuals/trizol\\_reagent.pdf](http://tools.invitrogen.com/content/sfs/manuals/trizol_reagent.pdf)