



# BLOCKING CLASS A SCAVENGER RECEPTORS

Created by: Chris Verschoor/updated by: Kyle Novakowski  
Updated: June 2012

Bowdish Lab, McMaster University  
Hamilton, ON, Canada  
[www.bowdish.ca](http://www.bowdish.ca)

## BACKGROUND

- This protocol details how to block class A scavenger receptors (CASR) in order to determine whether processes such as binding and phagocytosis is scavenger receptor mediated.

## EQUIPMENT

- Dextran sulphate, Sigma #D6001
- Chondroitin Sulfate, Sigma #C9819
- Fucoidan, Sigma #F563
- Polyinosinic acid, Sigma #P4154
- Polycytidylic acid, Sigma #P4903
- Macrophages
- Complete growth media (I.e. X-Vivo for monocyte derived macrophages)

## PROTOCOL

- 1) Prepare Class A Scavenger Receptor blockers and their respective negative controls at the following final concentrations:
  - o Dextran Sulfate (+) and Chondroitin Sulfate (-), 10 to 500  $\mu\text{g}/\text{ml}$  (100  $\mu\text{g}/\text{ml}$  optimal for monocyte-derived macrophages).
  - o Polyinosinic (+) and Polycytidylic acid (-), 10 to 500  $\mu\text{g}/\text{ml}$  (50  $\mu\text{g}/\text{ml}$  optimal for monocyte-derived macrophages).
  - o Fucoidan (no negative control), (-), 10 to 500  $\mu\text{g}/\text{ml}$  (100  $\mu\text{g}/\text{ml}$  optimal for monocyte-derived macrophages).
- 2) Incubate cells with blockers or negative controls for 30 mins at 37°C in complete media prior to stimulations or infections.
- 3) It is not necessary to wash off blockers prior to stimulation, however it may be warranted depending on the particular research question

## LINKS AND REFERENCES

- [www.bowdish.ca/lab/protocols](http://www.bowdish.ca/lab/protocols)