



SAMPLE PREPARATION FOR RNA-SEQ

Keith Lee & Michelle Pinto: April 2013

Bowdish Lab, McMaster University
Hamilton, ON, Canada

www.bowdish.ca

BACKGROUND

- For preparation of RNA samples for RNA-seq

NOTES

- Wear goggles, gloves and coat

EQUIPMENT

- Equipment:
 - o Vortex
 - o Centrifuge
 - o Pipettes and Tips, Tubes
 - o Aligent 2100 Bioanalyzer
- Materials:
 - o Aligent RNA 6000 Kit
 - o RiboZero Magnetic Kit – Human/Mouse/Rat (Epicentre MRZH11124)
 - o RNA Clean XP Kit
 - o Turbo DNase
 - o Superscript III First-Strand Synthesis System
 - o RNase H
 - o NEB Buffer 2
 - o Klenow fragment DNA polymerase
 - o Ampure XP
 - o 100% Ethanol
 - o RNase-free water

PROTOCOL

1. Check RNA integrity using the Aligent 2100 Bioanalyzer (get training from Jane-Ann)
 - a. Preparatory work
 - i. Aliquot 1 μ L of each RNA sample into a new 1.5mL microcentrifuge tube
 - ii. Sign-in on the sign-in sheet located beside the Bioanalyzer
 - iii. Remove reagents from fridge/freezer and thaw at room temperature for 30 minutes
 - iv. Turn the heat block on to 70°C
 - v. Clean Bioanalyzer electrodes
 1. Fill “RNase Zap” electrode cleaner with 350 μ L RNase Zap; shake chip

2. Open the machine lid and place electrode cleaner into the Bioanalyzer in the correct position; gently close lid using two hands and leave for 1 minute
 3. Fill "dH₂O" electrode cleaner with 350μL of dH₂O; shake chip
 4. Open the machine lid, remove the "RNase Zap" chip and insert the "dH₂O" chip; gently close lid using two hands and leave for 1 minute
 5. Open the machine lid, remove the "dH₂O" chip
 6. Keep lid open for 30 seconds for the water on the electrodes to evaporate before closing the lid
- b. Prepare RNA 6000 Nano Gel
 - i. Pipette 500μL of RNA 6000 nano gel matrix into spin filter
 - ii. Centrifuge gel matrix at 1500g for 10 minutes at room temperature
 - iii. Aliquot 65μL into a new 0.5mL tube
 - iv. The remainder can be stored at 4°C and can be used within 1 month of preparation
 - c. Prepare gel-dye mix
 - i. Vortex the RNA 6000 nano dye for 10 seconds
 - ii. Add 1μL of RNA 6000 nano dye to the 65μL filtered gel
 - iii. Vortex solution well and spin at 13000g for 10 minutes at room temperature
 - d. Prepare samples and ladder
 - i. Heat denature the samples (1μL) and ladder (1μL) at 70°C for 2 minutes
 - ii. Quick spin all samples to remove condensation on the lid
 - iii. Keep samples at room temperature until ready to use
 - e. Load the gel-dye mix
 - i. Place a new RNA 6000 nano chip on the chip priming station
 - ii. Pipette 9μL of the gel-dye mix into the well marked (G)
 - iii. Press down plunger until it is held by the clip after closing the priming station
 - iv. Wait 30 seconds and then release clip (plunger should move at least 0.3mL upwards within 1 second; if not, the chip may not be properly pressurized and may need to be reprimed)
 - v. Pipette 9μL of the gel-dye mix into the other 2 "G" wells
 - f. Load the RNA 6000 nano marker
 - i. Pipette 5μL of RNA 6000 nano marker into each sample well and the ladder well
 - g. Load the ladder and the samples
 - i. Add 1μL of ladder into the ladder well
 - ii. Add 1μL of sample into each sample well and 1μL of the marker into each unused well
 - iii. Vortex the loaded chip on IKA vortexer for 1 minute at 2400rpm
 - iv. Run the chip in the Aligent 2100 bioanalyzer within 5 minutes
 - v. Save results onto computer
 - h. Cleanup
 - i. Remove chip immediately after run is complete
 - ii. Fill "dH₂O" electrode cleaner with 350μL of dH₂O; shake chip
 - iii. Open the machine lid and insert the "dH₂O" chip; gently close lid using two hands and leave for 1 minute
 - iv. Open the machine lid, remove the "dH₂O" chip
 - v. Keep lid open for 30 seconds for the water on the electrodes to evaporate before closing the lid

- vi. Empty contents of the “RNase Zap” and “dH₂O” chips into the garbage but save the chips to reuse next time
 - vii. Turn off the Bioanalyzer, computer and heating block
 - viii. Complete the sign-in sheet with remaining information
 - i. Analysis
 - i. Ensure RNA has a high RNA Integrity Number (preferably >8)
2. rRNA depletion using RiboZero kit
- a. Bead-washing procedure
 - i. Pipette 225µL of magnetic beads (after vortexing) into a 1.5mL tube, per sample
 - ii. Place tube on magnetic stand for 1 minute
 - iii. Remove supernatant
 - iv. Add 225µL of RNase-free H₂O; vortex or pipette up and down
 - v. Place tube on magnetic stand for 1 minute
 - vi. Remove supernatant
 - vii. Remove tube from magnetic stand
 - viii. Add 65µL of magnetic bead resuspension solution and vortex
 - ix. Store at room temperature until ready
 - b. Treatment of RNA sample
 - i. Determine volume of RNA to use
 - 1. Based upon the highest RNA concentration in your samples, calculate a volume that contains either 1-2.5µg or 2.5-5µg of RNA
 - ii. In a RNase-free tube, mix:

x µL	RNase-free H ₂ O
4 µL	RiboZero reaction buffer
y µL	RNA sample (max 28µL if input is 1-2.5µg, max 26µL if input is 2.5-5µg)
z µL	RiboZero RNA removal sol'n (8µL for 1-2.5µg RNA, 10µL for 2.5-5µg RNA)
40 µL	Total
 - iii. Mix by pipetting gently
 - iv. Incubate at 68°C on a heat block for 10 minutes
 - v. Store at room temperature for 15 minutes
 - c. rRNA removal
 - i. Add 40µL (all) of treated RNA sample to magnetic beads and, without changing pipette tips, pipette up and down 10x, then vortex sample
 - ii. Place sample at room temperature for 10 minutes, mixing every couple of minutes
 - iii. Vortex at max speed for 5s
 - iv. Place at 50°C (on heating block) for 7 minutes (avoid significant condensation)
 - v. Place on magnetic stand for 1 minute
 - vi. Carefully remove supernatant (containing RNA, around 85-90µL) and transfer to new microcentrifuge tube
3. RNA Clean XP Kit purification
- a. Add 180µL of RNA clean beads to depleted RNA
 - b. Mix by pipetting up and down 10x gently
 - c. Incubate at room temperature for 15 minutes
 - i. During incubation, prepare fresh 80% EtOH
 - d. Place tube on magnetic stand for 5 minutes
 - e. Remove and discard supernatant, do not disturb the bead

- f. While still on stand, add 200 μ L of freshly prepared 80% EtOH
 - g. Incubate for 30s on stand, then remove supernatant without disturbing beads
 - h. Repeat 80% EtOH wash
 - i. Remove supernatant, allow beads to dry on stand for 10 minutes
 - j. Take tube off stand, add 32 μ L H₂O to elute RNA
 - k. Thoroughly resuspend beads by gently pipetting 10 times
 - l. Incubate at room temperature for 2 minutes
 - m. Place on magnetic stand for 5 minutes (or until solution is clear)
 - n. While still on stand, collect supernatant into a new tube and place on ice
 - o. Aliquot 2 μ L from each 32 μ L sample into a new tube for use with the Bioanalyzer (give samples to Christine in the Farncombe institute to analyze using the Pico kit)
 - p. Continue or store at -80°C
4. Turbo DNase digestion
 - a. Add the following to a RNase-free tube:

30 μ L	rRNA-depleted RNA (collected from previous step)
1.5 μ L	Turbo DNase
3.5 μ L	10x buffer
 - b. Mix contents with pipette
 - c. Incubate for 30 minutes at 37°C
 5. RNA Clean XP Kit purification (same as Step 3)
 - a. Same as step 2 but with 70 μ L beads (note: 2x volume of reaction)
 - b. Elute with 12 μ L H₂O
 - c. Continue or freeze at -80°C
 6. 1st strand cDNA synthesis
 - a. Add the following to a RNase-free PCR tube:

10 μ L	RNA from previous step
1 μ L	dNTPs
1 μ L	Random hexamers
12 μ L	Total
 - b. Place at 65°C on thermal cycler for 5 minutes
 - c. In order, prepare the following mix (per sample):

2 μ L	10x RT Buffer
4 μ L	MgCl ₂
2 μ L	DTT
1 μ L	RNase Out
1 μ L	SSC III
10 μ L	Total
 - d. Add the 10 μ L mix to the first mixture
 - e. Incubate for 10 minutes at room temperature
 - f. Incubate for 50 minutes at 50°C and then 5 minutes at 85°C in a thermal cycler
 - g. Continue or freeze at -80°C
 7. Ampure XP DNA purification
 - a. To the 22 μ L of the SSCIII mix from the previous step, add 44 μ L Ampure XP beads
 - b. Mix by pipetting 10 times
 - c. Incubate at room temperature for 5 minutes

- i. Prepare fresh 80% EtOH in the meantime
 - d. Place tubes on magnetic stand and let stand for 1 minute
 - e. Remove supernatant from each sample and discard
 - f. Add 200 μ L of 80% EtOH to each sample and incubate for 30s
 - g. Remove supernatant and repeat 200 μ L 80% EtOH wash
 - h. Air dry beads for 10-15 minutes (preferably until you can see cracks in the bead)
 - i. Remove tubes from the magnetic stand and elute DNA by adding 22 μ L H₂O to beads
 - j. Mix well by pipetting
 - k. Place tubes back on stand until solution is clear
 - l. Transfer 22 μ L supernatant to a new tube
 - m. Continue or freeze at -80°C
- 8. 2nd strand cDNA synthesis
 - a. Add the following to a nuclease-free PCR tube:

22 μ L	1 st strand cDNA from previous step
1 μ L	RNase H
3 μ L	dNTPs (10mM)
3 μ L	#2NEB Buffer
1 μ L	Klenow fragment (DNA polymerase)
<hr/>	
30 μ L	Total
 - b. Incubate at 16°C for 2 hours in a thermal cycler
 - c. Continue or freeze at -80°C
- 9. Ampure XP DNA purification
 - a. Same as step 7 but with 60 μ L beads
 - b. Elute with 50 μ L H₂O
 - c. Continue or freeze at -80°C
- 10. Covaris sonication (Covaris S220)
 - a. Use the following protocol for 50 μ L samples, 150 base pairs
 - i. Peak power: 175W
 - ii. Duty factor: 20%
 - iii. Cycles/burst: 200
 - iv. Time: 155s
 - b. Continue or freeze at -80°C
- 11. End repair
 - a. Make master mix (per sample)

10 μ L	10x Buffer
5 μ L	NEB End Repair Enzyme Mix
35 μ L	Sterile (nuclease-free) H ₂ O
<hr/>	
50 μ L	Total
 - b. Add the 50 μ L master mix to each 50 μ L fragmented cDNA sample
 - c. Place in thermal cycler for 30 minutes at 20°C
- 12. Ampure XP DNA Purification
 - a. Same as step 7 but with 160 μ L beads (1.6x volume)
 - b. Elute with 42 μ L H₂O
 - c. Continue or freeze at -80°C

13. dA-tailing

- a. Add the following to a nuclease-free PCR tube:

42 μ L End-repaired cDNA
5 μ L dA-tailing reaction buffer
3 μ L Klenow fragment

50 μ L Total

- b. Place in thermal cycler for 30 minutes at 37°C

14. Ampure XP DNA purification

- a. Same as step 7 but with 90 μ L beads (1.8x volume)
b. Elute with 25 μ L H₂O
c. Continue or freeze at -80°C

15. Adapter ligation

- a. Add the following to a nuclease-free PCR tube:

25 μ L dA-tailed cDNA
10 μ L Adapters for Illumina
10 μ L Ligase buffer (5x)
5 μ L Ligase

50 μ L Total

- b. Place in thermal cycler for 15 minutes at 20°C
c. Add 3 μ L of USER enzyme to each sample
d. Place in thermal cycler for 15 minutes at 37°C

16. Ampure XP DNA purification

- a. Same as step 7 but with 50 μ L beads (1.0x volume)
b. Elute with 34 μ L H₂O
c. Continue or freeze at -80°C

17. PCR amplification with indexes

- a. Create the following master-mix (per sample)

0.25 μ L Taq polymerase
5 μ L 10x Taq buffer
1.5 μ L MgCl₂ (50mM)
1 μ L dNTPs (10mM)
2.5 μ L Universal Primer (25 μ M)
3.25 μ L H₂O

13.5 μ L Total

- b. Add the master mix to the 34 μ L cDNA sample

- c. Add 2.5 μ L of a unique index primer (25 μ M) to each sample

- i. For the index primers, use NEBNext Index Sets 1 and 2, or design custom primers

- d. Amplify samples with the following PCR reaction:

i. Initial denaturation 94°C 3 minutes
ii. Denaturation 94°C 45 seconds
iii. Annealing 65°C 30 seconds
iv. Extension 72°C 30 seconds
v. Final extension 72°C 10 minutes
vi. Hold 4°C

} 4-8 cycles

18. Ampure XP DNA purification

- a. Same as step 7 but with 50 μ L beads (1.0x volume)
 - b. Elute with 27 μ L H₂O
 - c. Continue or freeze at -80°C
19. Quality Check and RNA sequencing run
- a. Provide samples to Christine King in the Farncombe Metagenomics Facility, who will:
 - i. Run a Bioanalyzer run (DNA chip) to validate the size of your cDNA library (fragments are expected to be approximately ~260 base pairs) and to rule out the possibility of adapter-dimers in your library
 - ii. Run a qPCR reaction to quantify the concentration of your samples
 - iii. Sequence your samples using an Illumina HiSeq machine