



FLOW CYTOMETRY STAINING

Charles Yin

Last Updated: 21 August 2014

Bowdish Lab, McMaster University

Hamilton, ON, Canada

www.bowdish.ca

EQUIPMENT AND MATERIALS

- Equipment:
 - o Centrifuge set at 4°C
 - o Bucket with ice
 - o PBS
- Materials:
 - o Antibodies for flow cytometry
 - o FACS wash (0.5% BSA, 5mM EDTA in PBS)
 - o FACS block (1% BSA, 5% FBS in FACS wash)
 - o 4% paraformaldehyde (PFA) in PBS

PROTOCOL

1. Wash cells with cold PBS.
2. Lift cells and re-suspend in 1ml cold PBS.
 - a. Cells should not be exposed to trypsin nor scraped off the plate. Use cell lifters where possible.
 - b. Keep cells cold on ice at all points after this step to prevent membrane turnover.
3. Pellet cells by centrifuging at 4°C at 2,000rpm for 10min.
4. Re-suspend pellet in 100µl FACS block.
5. Incubate on ice for 1h.
6. Add primary and secondary (or conjugated) antibodies.
 - a. Most protocols recommend adding primary and secondary antibodies separately. I have seen no difference either way. Adding both antibodies together reduces the total time of this protocol.
 - b. Keep cells away from light as much as possible after this step to avoid photobleaching of fluorophores attached to the antibodies.
7. Incubate on ice for 1h.
8. Wash 3× with 1ml cold PBS.
9. If cells are to be analyzed immediately, re-suspend cells in 200ul cold PBS. Otherwise, re-suspend cells in 200ul of cold 4% PFA and store up to 24h at 4°C away from light.
10. Strain cells into flow cytometry tubes prior to analysis. Keep tubes covered when not in use.