

Scavenger Receptor Genotyping Protocol

Mike Dorrington - Feb 12, 2011

Materials:

1. Mouse Tails
2. Extract 'n' Amp Tissue PCR Kit (Sigma)
3. Primer mix (10 μ M each)
4. Agarose gel

Procedure:

1. For each reaction (each tail) mix together 100 μ L of extraction solution and 25 μ L of tissue prep solution (from Extract n Amp kit)
2. In a 15mL eppendorf tube, add 125 μ L of mix to each tail end with the cut end facing down
 1. Room temperature for 10min
 2. Then 95 $^{\circ}$ C water bath for 3 min
3. Add 100 μ L of Neutralization solution B and vortex liberally
4. The tails (in solution) can either be used for PCR now or be stored at 4 $^{\circ}$ C until you want to perform the PCR reaction (can be stored this way for up to 6 months)
5. For each PCR reaction mix together:
 1. 4.5 μ L PCR-grade H₂O
 2. 10 μ L PCR mix (from kit)
 3. 1.5 μ L primer mix (forward and back)
 4. 4 μ L of tissue extract solution
6. Mix gently and perform PCR cycling in a thermal cycler:
 1. Denature 3min at 94 $^{\circ}$ C
 2. Then 30 cycles of:
 1. Denature 1min at 94 $^{\circ}$ C
 2. Anneal 1min at 56 $^{\circ}$ C
 3. Extend 1min at 72 $^{\circ}$ C
 3. Final extension 10min at 72 $^{\circ}$ C
 4. Hold at 4 $^{\circ}$ C indefinitely
7. The PCR product can now be directly loaded into an agarose gel and run as per usual

Primer Sets:

SRA WT (~250bp)

F: 5' ACC TTA TAG ACA CGG GAC GCT TCC AGA A 3'

R: 5' GAC TCT GAC ATG CAG TGT TTC TGT A 3'

SRA KO (~350bp)

F: 5' ACC TTA TAG ACA CGG GAC GCT TCC AGA A 3'

R: 5' AGG AGT AGA AGG TGG CGC GAA GG 3'

MARCO WT (~500bp)

F: 5' CAG CTG GGT CCA TAC CAG C 3'

R: 5' CTG GAG AGC CTC GTT CAC C 3'

MARCO KO (~850bp)

F: 5' CCA CGC TCA TCG ATA ATT TCA C 3'

R: 5' GCC TGC AGT GGC CGT CGT TTT A 3'