



COOMASSIE BLUE STAIN

Created by: Jason Fan
April 9, 2016

Bowdish Lab, McMaster University
Hamilton, ON, Canada
www.bowdish.ca

Preparation:

1. Remove 1.5M Tris (pH 8.8) (w/ SDS) from 4°C fridge
2. Remove 1.0M Tris (pH 6.8) (w/ SDS) from 4°C fridge
3. Remove 30% acrylamide mix from 4°C fridge
4. TEMED
5. dH₂O
6. 70% ethanol
7. Coomassie dye recipe (the order of preparation is critical):
 - a. Dissolve 25g of Aluminum Sulfate in 200 mL of Millipore water to create a 5% w/v concentration
 - b. Add 50 mL of 100% ethanol for a final concentration of 10% v/v
 - c. Dissolve 0.1 g of Coomassie Brilliant Blue G-250 (Sigma) to create a 0.02% w/v concentration; immediately mix well by swirling and inverting the bottle repeatedly
 - d. Add 12 mL of concentrated o-phosphoric acid for a final concentration of 2% v/v. The addition of the acid to the ethanol will facilitate colloid production.
 - e. Fill up to 500 mL with Millipore water

Protocol

1. Assemble gel cassette
 - a. Place gray rubber strip on the bottom of one of the two 'clip on' stations (large plastic 'stations')
 - b. Take two glass plates (one has to be 1.0mm plate) and wash them gently with soap water
 - i. Wear gloves
 - ii. Hold plate by the edges
 - iii. Allow to air dry on a paper towel
 - c. Place the two plates together and insert them into the green holder
 - i. Secure by snapping the sides
 - ii. Ensure the bottom is even to prevent leaking
 - d. 'clip' the 1.0mm glass plate/green holder unit onto the plastic station
 - e. Add 5 mL of water to test for leakage
 - i. Dump water
 - ii. Dry gently with Kim wipe
2. Prepare Separating Gel (make 1.5x the amount to be safe)
 - a. 1.5 mL dH₂O
 - b. 2.3 mL 30% acrylamide

- c. 1.3 mL 1.5M Tris (pH 8.8) (w/ SDS)
 - d. 40 μ l 10% APS
 - i. Weigh 0.1 g of APS powder and dissolved in 1 mL of dH₂O
 - ii. Shake solution gently
 - e. 2 μ l TEMED
 - i. Shake solution
3. Add gel to cassette
 - a. Add ~5 mL of Separating Gel to cassette
 - b. Draw up remaining in a 1 mL pipet
 - c. Wait 20 minutes (could take longer) for polymerization to occur
 - i. Test by regularly ejecting pipet to see if gel has polymerized
 - d. Add 70% ethanol to cassette (up to the brim of the two glass plates)
4. After polymerization has occurred, pour out ethanol and wash with dH₂O. Dry with Kim wipe (gently) as best as possible
5. Prepare Stacking Gel (make 1.5x amount to be safe)
 - a. 863 μ l H₂O
 - b. 213 μ l 30% acrylamide mix
 - c. 156 μ l 1.0M Tris (pH 6.8) (w/ SDS)
 - d. 13 μ l 10% APS
 - i. Shake
 - e. 1 μ l TEMED
 - i. Shake
6. Add Stacking Gel to cassette
 - a. Add enough gel to reach brim of glass plates
 - i. Takes ~3 mL
 - b. Insert green comb straight into gel
 - c. ****turning on heat block now will save time – turn temperature to 95°C****
 - d. Gel will polymerize quickly – wait 5 min. (could take longer)
 - e. Once gel has polymerized, remove comb by lifting **vertically** and gently
7. Prepare samples
 - a. Prepare a DTT + 3x sample buffer solution
 - i. Add 250 μ l DTT (found in -20°C freezer in box labelled ‘common reagents’) **per** 750 μ l 3x SDS sample buffer (dark purple solution)
 - b. Add 50 μ l of DTT solution **per** 100 μ l of cell sample
 - c. Place sample on 95°C heat block for 10 minutes
8. Loading Samples
 - a. Make 1x running buffer
 - i. Add 240 mL of 5x stock solution to 960 mL MilliQ water for total volume of 1.2 L
 - b. Remove two glass plates with gel from green holder

- c. Insert glass plates into plastic holder that can hold two sets of plates
 - i. Pull apart the side locking handles
 - ii. Place glass plate into one of the sides (1.0 mm plate facing inward)
 - iii. Place a purple plastic plate on the other side (words facing inwards)
 - iv. Lock in place
 - d. Place plastic holder in tank with corresponding currents correct
 - e. Add 20 μ l of protein ladder (found in -20°C fridge in common reagents box) to first well
 - f. Add 30 μ l of each sample into subsequent wells
9. Running Samples
- a. Pour 1.2L of 1x running buffer into tank
 - b. Push bubbles away from the plastic holder if there are any
 - c. Place lid on with correct currents
 - d. Run sample for ~40 minutes at 200V
 - i. Check to see if small bubbles rise from bottom of the tank
10. Remove gel gently and place onto large Petri dish
- a. Wash twice with Millipore water to remove excess SDS
 - b. Submerge with required stain and place on shaker overnight
 - c. Bands should appear within 2 hours but overnight exposure will maximize the appearance of bands; the gel will have appeared to have shrunken but this is normal
11. Destain the gel by washing repeatedly with Millipore water