

# PHENOL-SULFURIC ACID MANNOSE DETECTION ASSAY

Jason Fan

November 9, 2015

Bowdish Lab, McMaster University

Hamilton, ON, Canada

[www.bowdish.ca](http://www.bowdish.ca)

## BACKGROUND AND NOTES

- This assay is sensitive to mannose from 1 – 100  $\mu\text{g}$  but is sensitive towards other carbohydrates
- This assay was used for the quantification of complex carbohydrates that constitutes bacterial cell wall
  - o Whole Gram-positive bacterial lysates can be detected using this assay providing a fast, colorimetric alternative to bacterial quantification
- Volumes suggested below can be manipulated but the **ratio of reagents must be maintained**

## EQUIPMENT


- Equipment:
  - o Spectrophotometer capable of measuring absorbance at 490nm
  - o 96-well plate
- Materials:
  - o Concentrated  $\text{H}_2\text{SO}_4$
  - o 5% (w/v) Phenol

## PROTOCOL

- Preparatory Work:
    - o Create a 5% w/v phenol in water solution; small volumes of phenol are required per sample but a stock solution may be made for convenience
      - Phenol is primarily non-polar but at low concentrations, it is soluble in water
    - o The following steps can be done right in a 96-well plate or in Eppendorf tubes; absorbance measurements are still to be done in a plate reader
1. Add 50  $\mu\text{l}$  of sample into 96 well plate
    - a. It is best to create triplicates per measured sample
    - b. Ensure the sample vehicle does not contain sugars
  2. **Forcible and quickly** eject 150 $\mu\text{l}$  concentrated  $\text{H}_2\text{SO}_4$  into each well
    - a. **This reaction is highly exothermic; be careful of fumes and heat**
    - b. Quickly ejecting into each well helps with mixing
  3. Gently shake the plate until solution looks homogenized
    - a. Sugars present in the sample will react with the acid to quickly change into a yellow-brown colour
  4. Add 30 $\mu\text{l}$  of 5% phenol
  5. Let stand for 10 minutes at room temperature

6. Read absorbance at a wavelength of 490 nm

## RESULTS



$$Y = 0.03128 * X + 0.01288$$

$$R^2 = 0.9981$$

**Figure #1: Phenol-Sulfuric Acid Mannose Detection Assay Standard Curve.** The standard curve displays weight of mannose detected. All values can be considered in concentrations by accounting the 50  $\mu$ l sample volume.

## LINKS AND REFERENCES

1. Masuko T, Minami A, Iwasaki N, Majima T, Nishimura S-I, Lee YC. Carbohydrate analysis by a phenol-sulfuric acid method in microplate format. *Analytical Biochemistry*. 2005 Apr;339(1):69-72.