

# WHOLE BLOOD MONOCYTE STAIN-FLOW CYTOMETRY

Created by: **Grace Teskey**      Date: **March 2017**  
Bowdish Lab, McMaster University  
Hamilton, ON, Canada  
[www.bowdish.ca](http://www.bowdish.ca)

## EQUIPMENT/MATERIALS

- eBioscience 1-Step fix/lyse solution (10X)
- Antibodies (See page 4)
- FACs Wash (0.5% (w/v) BSA, 5mM EDTA (pH 7.4-7.6), for 500mL 2.5g BSA, 5mL of 0.5M EDTA)
- 2mL microcentrifuge tubes
- 5mL polystyrene round bottom tubes
- 5mL polystyrene round bottom tubes with cell-strainer cap (Corning REF:352235)
- eBioscience OneComp eBeads
- CountBright™ Absolute Counting Beads
- FlowJo software

## PROTOCOL

### Cell Staining

1. Prepare lymphocyte stain in a 2mL microcentrifuge tube. Stain for 1 sample is made up to 50  $\mu$ L in FACs Wash (8.5  $\mu$ L and antibody + 41.5  $\mu$ L FACs Wash) and is in 3X working concentration. *Concentrations shown below are for 50 $\mu$ L final volume.*


Cell Surface Marker	Fluorophore	Antibody Concentration	Volume ( $\mu\text{L}$ )
CD45	BV510	1/100	0.5
CD16	PE-Cy7	1/100	0.5
CD14	BV421	1/100	0.5
CCR2	PE	1/50	1
CD11b	APC	1/50	1
HLA-DR	PerCPCy5.5	1/100	0.5
CD15	BV650	1/100	0.5
CX3CR1	FITC	1/50	1
CD3	AF700	1/50	1
CD56	AF700	1/50	1
CD19	AF700	1/50	1

2. Add 100  $\mu\text{L}$  of whole blood that was collected in a heparinized or EDTA coated blood collection tube. Incubate for 30 minutes at room temperature away from light.
3. Top up 2mL microcentrifuge tube containing stain and blood with 1-Step fix/lyse solution (made to 1X in distilled water). Incubate for 10 minutes at room temp away from light with intermittent inversion. *This step is for RBC lysis as well as fixation of white blood cells.*
4. Spin down at 2000rpm for 5 minutes at 4°C
5. Aspirate supernatant. Wash pellet with 2mL FACs Wash. Repeat step 4.
6. Aspirate supernatant and resuspend in 240  $\mu\text{L}$  of FACs Wash. *If running the sample on cytometer immediately, continue to step 7. If not, store samples at 4 °C away from light for up to 24 hours)*
7. Filter cells through 0.45  $\mu\text{m}$  mesh to ensure single cell suspension and add 10  $\mu\text{L}$  of count beads to each sample if counts are needed.

### Compensation Controls- To be made right before running flow

1. In a 5mL polystyrene tube, dilute 2 drops of OneComp eBeads to 2.5mL with FACs Wash and vortex. *Be sure to vortex OneComp eBeads before use.*
2. Aliquot 240 $\mu$ L of diluted eBeads into 10 polystyrene tubes.
3. Add 0.5 $\mu$ L of each antibody used in the stain into one tube (i.e. 0.5 $\mu$ L of CD3 antibody into tube 1, 0.5 $\mu$ L of CD4 antibody into tube 2, etc) leaving one tube unstained. Vortex tubes. *Only one compensation tube needed for AF700.*
4. Keep tubes away from light until needed.

### Monocyte Gating Strategy



**Antibody List**

Antibody	Fluorophore	Company	Clone	Catalog#
CD45	BV510	BioLegend	HI30	304036
CD16	PE-Cy7	eBioscience	eBioCB16	25-0168-42
CD14	BV421	BioLegend	M5E2	301830
CD11b	APC	BD Pharmingen	ICRF44	561015
CCR2	PE	BioLegend	K036C2	357206
HLADR	PerCPy5.5	eBioscience	LN3	45-9956-42
CD15	BV650	BioLegend	W6D3	323034
CX3CR1	FITC	BioLegend	2A9-1	D070-4
CD3	AF700	BD Biosciences	UCHT1	557943
CD19	AF700	eBioscience	HIB19	56-0199-42
CD56	AF700	BioLegend	HCD56	318316