



# WHOLE BLOOD MONOCYTE STAIN-FLOW CYTOMETRY

Created by: Grace Teskey

Date: March 2017

Updated: July 2018

Bowdish Lab, McMaster University  
Hamilton, ON, Canada

[www.bowdish.ca](http://www.bowdish.ca)

## EQUIPMENT/MATERIALS

- eBioscience 1-Step fix/lyse solution (10X) (00-5333-57)
- Antibodies (See page 3)
- FACs Wash (0.5% (w/v) BSA, 5mM EDTA (pH 7.4-7.6), for 500mL 2.5g BSA, 5mL of 0.5M EDTA)
- 2mL microcentrifuge tubes
- 5mL polystyrene round bottom tubes
- 5mL polystyrene round bottom tubes with cell-strainer cap (Corning REF:352235)
- eBioscience OneComp eBeads (01-1111-41)
- CountBright™ Absolute Counting Beads (C36950)
- FlowJo software

## PROTOCOL

### Buffer Preparation

1. Prepare a 1x working solution of 1-Step Fix/Lyse buffer by mixing 1 part of 10x Fix/Lyse solution with 9 parts of distilled water.

### Cell Staining

1. Prepare lymphocyte stain in a 2mL microcentrifuge tube. Stain for 1 sample is made up to 50  $\mu$ L in FACs Wash (8.25  $\mu$ L and antibody + 41.75  $\mu$ L FACs Wash) and is in 3X working concentration. *Concentrations shown below are for 50 $\mu$ L final volume.*

Cell Surface Marker	Fluorophore	Antibody Concentration	Volume ( $\mu$ L)
CD45	BV510	1/100	0.5
CD16	PE-Cy7	1/100	0.5
CD14	BV421	1/100	0.5
CCR2	PE	1/50	1.0
CD11b	APC	1/50	1.0
HLA-DR	PerCPCy5.5	1/100	0.5
CD15	BV650	1/200	0.25
CX3CR1	FITC	1/50	1.0
CD3	AF700	1/50	1.0
CD56	AF700	1/50	1.0
CD19	AF700	1/50	1.0

**Note:** for each sample, four stains are required. One unstained control (50 $\mu$ L of FACS wash). One CD16 FMO stain (all antibodies in the cocktail except for CD16-PE-Cy7). One CCR2 ISotype control (all antibodies in the cocktail and replace CCR2-PE with an ISotype antibody) and one complete stain as described above.

2. Add 100  $\mu$ L of whole blood that was collected in a heparinized or EDTA coated blood collection tube. Incubate for 30 minutes at room temperature away from light.
3. Top up 2mL microcentrifuge tube containing stain and blood with 1-Step fix/lyse solution (made to 1X in distilled water). Incubate for 10 minutes at room temp away from light with intermittent inversion. *This step is for RBC lysis as well as fixation of white blood cells.*
4. Spin down at 2000rpm for 5 minutes at 4°C.
5. Aspirate supernatant. Wash pellet with 2mL FACS Wash. Spin down at 2000rpm for 5 minutes at 4°C.

6. Aspirate supernatant and resuspend in 240  $\mu$ L of FACs Wash. *If running the sample on cytometer immediately, continue to step 7. If not, store samples at 4  $^{\circ}$ C away from light for up to 24 hours)*
7. Filter cells through 0.45  $\mu$ m mesh to ensure single cell suspension and add 10  $\mu$ L of count beads to each sample if counts are needed.

**Compensation Controls-** *To be made right before running flow*

1. In a 5mL polystyrene tube, dilute 2 drops of OneComp eBeads to 2.5mL with FACs Wash and vortex. *Be sure to vortex OneComp eBeads before use.*
2. Aliquot 240 $\mu$ L of diluted eBeads into 10 polystyrene tubes.
3. Add 0.5 $\mu$ L of each antibody used in the stain into one tube (i.e. 0.5 $\mu$ L of CD45 antibody into tube 1, 0.5 $\mu$ L of CD16 antibody into tube 2, etc) leaving one tube unstained. Vortex tubes. *Only one compensation tube needed for AF700.*
4. Keep tubes away from light until needed.

**Antibody List**

Antibody	Fluorophore	Company	Clone	Catalog#
CD45	BV510	BioLegend	HI30	304036
CD16	PE-Cy7	eBioscience	eBioCB16	25-0168-42
CD14	BV421	BioLegend	M5E2	301830
CD11b	APC	BD <sup>Pharmingen</sup>	ICRF44	561015
CCR2	PE	BioLegend	K036C2	357206
HLADR	PerCPy5.5	eBioscience	LN3	45-9956-42
CD15	BV650	BioLegend	W6D3	323034
CX3CR1	FITC	BioLegend	2A9-1	D070-4
CD3	AF700	BD <sup>Bioscience</sup>	UCHT1	557943
CD19	AF700	eBioscience	HIB19	56-0199-42
CD56	AF700	BioLegend	HCD56	318316

# Monocyte Gating Strategy

Monocyte Gating Strategy Updated: July 2018

