

WHOLE BLOOD T CELL SENESENCE-FLOW CYTOMETRY

Created by: Allison Kennedy Date: February 2019
Bowdish Lab, McMaster University
Hamilton, ON, Canada
www.bowdish.ca

EQUIPMENT/MATERIALS

- eBioscience 1-Step fix/lyse solution (10X)
- Antibodies (See page 5)
- FACs Wash (0.5% (w/v) BSA, 5mM EDTA (pH 7.4-7.6), for 500mL 2.5g BSA, 5mL of 0.5M EDTA)
- 2mL microcentrifuge tubes
- 5mL polystyrene round bottom tubes
- 5mL polystyrene round bottom tubes with cell-strainer cap (Corning REF:352235)
- eBioscience OneComp eBeads
- CountBright™ Absolute Counting Beads
- FlowJo software

PROTOCOL

Cell Staining

1. Prepare lymphocyte stain in a 2mL microcentrifuge tube. Stain for 1 sample is made up to 50 μ L in FACs Wash (9.25 μ L and antibody + 40.75 μ L FACs Wash) and is in 3X working concentration. *Concentrations shown below are for 50 μ L final volume.*

| Cell Surface Marker | Fluorophore | Antibody Concentration | Volume (μ L) |
|---------------------|-----------------|------------------------|-------------------------|
| CD19 | eFluor 450 | 1/50 | 1.0 |
| CD3 | PerCPCy5.5 | 1/200 | 0.25 |
| CD4 | APC.Cy7 | 1/200 | 0.25 |
| CD8 | BV510 | 1/200 | 0.25 |
| CD45 | Alexa Fluor 700 | 1/100 | 0.5 |
| CD45RA | FITC | 1/100 | 0.5 |
| CD14 | eFluor 450 | 1/75 | 1.5 |
| CD15 | eFluor 450 | 1/100 | 0.5 |
| CCR7 | PE-eFluor 610 | 1/50 | 1 |
| CD57 | PE | 1/50 | 1 |
| CD28 | APC | 1/50 | 1 |
| CD56 | eFluor 450 | 1/100 | 0.5 |
| Zombie Yellow | Qdot585 | 1/50 | 1 |
| | | | Ab=9.25 FACS = 40.75 |


2. Add 100 μ L of whole blood that was collected in a heparinized or EDTA coated blood collection tube. Incubate for 30 minutes at room temperature away from light.
3. Top up 2mL microcentrifuge tube containing stain and blood with 1-Step fix/lyse solution (made to 1X in distilled water). Incubate for 10 minutes at room temp away from light with intermittent inversion. *This step is for RBC lysis as well as fixation of white blood cells.*
4. Spin down at 2000rpm for 5 minutes at 4°C
5. Aspirate supernatant. Wash pellet with 2mL FACs Wash. Repeat step 4.

6. Aspirate supernatant and resuspend in 240 μL of FACs Wash. *If running the sample on cytometer immediately, continue to step 7. If not, store samples at 4 $^{\circ}\text{C}$ away from light for up to 24 hours)*
7. Filter cells through 0.45 μm mesh to ensure single cell suspension and add 10 μL of count beads to each sample if counts are needed.

Compensation Controls- *To be made right before running flow*

1. In a 5mL polystyrene tube, dilute 2 drops of OneComp eBeads to 2.5mL with FACs Wash and vortex. *Be sure to vortex OneComp eBeads before use.*
2. Aliquot 240 μL of diluted eBeads into 10 polystyrene tubes.
3. Add 0.5 μL of each antibody used in the stain into one tube (i.e. 0.5 μL of CD3 antibody into tube 1, 0.5 μL of CD4 antibody into tube 2, etc) leaving one tube unstained. Vortex tubes. *Only one compensation tube needed for eFluor 450.*
4. Keep tubes away from light until needed.

Senescent T cell Gating (Below example is with Isolated PBMCs):



Antibody List

| Marker | Fluorochrome | Concentration | Clone | Manufacturer | Ref/Cat # |
|---------------|-----------------|---------------|--------|--------------|------------|
| CD19 | eFluor 450 | 1/50 | HIB19 | eBioscience | 48-0199-41 |
| CD3 | PerCPCy5.5 | 1/200 | OKT3 | eBioscience | 45-0037-41 |
| CD4 | APC.Cy7 | 1/200 | OKT4 | BioLegend | 317417 |
| CD8 | BV510 | 1/200 | SK1 | BioLegend | 344731 |
| CD45 | Alexa Fluor 700 | 1/100 | 2D1 | eBioscience | 56-9459-41 |
| CD45RA | FITC | 1/100 | JS-83 | eBioscience | 11-9979-41 |
| CD14 | eFluor 450 | 1/75 | 61D3 | eBioscience | 48-0149-41 |
| CD15 | eFluor 450 | 1/100 | MMA | eBioscience | 48-0158-41 |
| CCR7 | PE-eFluor 610 | 1/50 | 3D12 | eBioscience | 61-1979-41 |
| CD57 | PE | 1/50 | TB01 | Invitrogen | 12-0577-41 |
| CD28 | APC | 1/50 | CD28.2 | eBioscience | 17-0289-41 |
| CD56 | eFluor 450 | 1/100 | TULY56 | eBioscience | 48-0566-41 |
| Zombie Yellow | Qdot585 | 1/50 | -- | ThermoFisher | L34967 |