



TNF GENOTYPING

Created/updated by: Jessica Wallace Date: August 7, 2019
Bowdish Lab, McMaster University
Hamilton, ON, Canada
www.bowdish.ca

BACKGROUND

This protocol is used to genotype TNF knockout ($TNF^{-/-}$), TNF heterozygote ($TNF^{+/-}$), and wild type ($TNF^{+/+}$) mice from the Bowdish Lab Litter Mate Heterozygote (LMH) breeding colony.

NOTES

- Expected results
 - Wild type: 183 bp
 - TNF mutant: 318 bp
 - Heterozygote: 183 bp and 318 bp
- Protocol adapted from Stock Number: 003008 from the Jackson Laboratory
https://www2.jax.org/protocolsdb/f?p=116:5:0::NO:5:P5_MASTER_PROTOCOL_ID,P5_JRS_CODE:22433,03008
- Additional PCR tips can be consulted at: <https://www.jax.org/jax-mice-and-services/customer-support/technical-support/genotyping-resources/general-pcr-reminders-dos-and-donts>

MATERIALS

- PCR grade water
- 10X PCR buffer (Invitrogen, 10342020)
- 50 μ M MgCl₂ (Invitrogen, 10342020)
- 10 μ M dNTP mix (ThermoFisher Scientific, 18-427-088)
- 100 μ M Common TNF genotyping forward primer (IDT, oIMR4182, 5' TAG CCAGGA GGG AGA ACA GA 3')
- 100 μ M Wild type reverse TNF genotyping primer (IDT, oIMR4183, 5' AGT GCC TCT TCT GCC AGT TC 3')
- 100 μ M Mutant reverse TNF genotyping primer (IDT, oIMR7297, 5' CGT TGG CTA CCC GTG ATA TT 3')
- Genomic DNA
- Taq DNA polymerase, recombinant (ThermoFisher Scientific, 10342020)
- 8-Strip PCR Tube (DiaMed, DIATEC420-1377)
- 100bp + 3K DNA ladder (SmoBio, DM2300)
- 6x Purple gel loading dye (New England Biolabs, B7024S)

PROTOCOL

1. Prepare a PCR master mix comprised of PCR grade water, 10X PCR buffer, 50 μ M MgCl₂, 10 μ M dNTPs, 10 μ M common TNF genotyping forward primer, 10 μ M wild type reverse TNF genotyping primer, and 10 μ M mutant reverse TNF genotyping primer. Add Taq polymerase last. Mix by vortexing. Keep on ice.
 - a. A single PCR reaction contains the following components at the following volumes:
 - i. PCR grade water 16.3 μ l

ii.	10X PCR buffer	2.5µl
iii.	50µM MgCl ₂	0.75µl
iv.	10µM dNTP mix	0.50µl
v.	Common TNF forward genotyping primer	1.25µl
vi.	Wild type reverse TNF genotyping primer	1.25µl
vii.	Mutant reverse TNF genotyping primer	1.25µl
viii.	Genomic DNA	1µl
ix.	Taq DNA polymerase	0.2µl

2. Add 24µl of the PCR master mix to sterile 8-strip 0.2mL PCR tubes
 - a. Mix by vortexing with each new 8-strip
 - b. Keep 8-strip 0.2mL PCR tubes on ice

3. Add 1µl genomic DNA to sterile 8-strip 0.2mL PCR tubes containing PCR master mix
 - a. Gently mix by flicking the 8-strip 0.2mL PCR tubes
 - b. Keep 8-strip 0.2mL PCR tubes on ice

4. Briefly pulse centrifuge 8-strip 0.2mL PCR tubes containing genomic DNA and PCR master mix

5. Place samples in thermocycler (Eppendorf, Mastercycler, Pro S) and run the following touchdown amplification program:
 - a. 94°C 5 minutes
 - b. 94°C 45 seconds
 - c. 65°C 30 seconds
 - d. 68°C 40 seconds x10 cycles: T_m decreases by 0.5° every cycle
 - e. 94°C 45 seconds
 - f. 55°C 30 seconds
 - g. 72°C 40 seconds x35 cycles
 - h. 72°C 7 minutes
 - i. 4°C Infinite hold

6. Remove samples from thermocycler
 - a. Add 4µl 6x purple gel loading dye to each sample
 - b. Gently mix by flicking the 8-strip 0.2mL PCR tubes
 - c. Keep 8-strip 0.2mL PCR tubes on ice

7. Briefly pulse centrifuge 8-strip 0.2mL PCR tubes containing PCR products and loading dye

8. Separate PCR products by gel electrophoresis on a 1.5% agarose gel
 - a. Small gel:
 - i. 50mL 1X TAE buffer
 - ii. 0.75g agarose
 - iii. 2.5µl RedSafe nucleic acid staining solution
 - b. Medium gel:
 - i. 100mL 1X TAE buffer
 - ii. 1.5g agarose

- iii. 5 μ l RedSafe nucleic acid staining solution
- c. Run at: Voltage: 120V, Duration: 25 minutes
- d. Large gel:
 - i. 350mL 1X TAE buffer
 - ii. 5.25g agarose
 - iii. 15 μ l RedSafe nucleic acid staining solution
- e. Run at: Voltage: 120V, Duration: 60-75 minutes