



MONOCYTE CELL SORTING FOR RNASEQ

Created by: Allison Kennedy Date: June 2019

Bowdish Lab, McMaster University
Hamilton, ON, Canada
www.bowdish.ca

EQUIPMENT/MATERIALS

- Antibodies (See page 4)
- FACs Wash (0.5% (w/v) BSA, 5mM EDTA (pH 7.4-7.6))
- eBioscience RBC Lysis Buffer (10X) (420301)
- RPMI Media
- PBS
- Leucosep tubes filled with ficoll (VWR: 89048-934)
- 50 mL conical tubes
- 5mL polystyrene round bottom tubes (FACS Tubes)
- 5mL polystyrene round bottom tubes with cell-strainer cap (352235)
- 5mL polypropylene round bottom tubes (FACS Tubes)
- BD Pharmingen 7 AAD (559925)
- eBioscience OneComp eBeads (01-1111-41)
- Qiagen QIAzol Lysis Reagent (79306)

PROTOCOL

Cell Sorting for Human Monocytes Comparing PBMCs to Fresh Blood for RNA seq.

*** from 100uL of whole blood average of 8000 classical, 230 Intermediate, 300 Non-Classical Monocytes*

*** from 30 million PBMCs average 6.29×10^5 classical, 0.24×10^5 intermediate and 0.22×10^5 non classical*

**** aiming for 50 000 x all three classifications.*

WHOLE BLOOD MONOCYTES

- 1) Collect 3x10 mL of whole blood in a heparin collection tube. Allow to sit at room temperature for at least 30 minutes before processing.
- 2) Spin at 1500 rpm for 10 mins and remove plasma.
- 3) Combine the volume of blood 1:1 with PBS in a 50mL conical tube.
- 4) Divide the volume into two Leucosep tubes and centrifuge without a break at 2300 rpm for 10 mins.
- 5) Remove plasma layer from the Leucosep tube and combine buffy layers from both tubes into one 50mL conical tube.
- 6) Wash cells by adding cold PBS up to 50mL. Centrifuge at 1500rpm for 5 minutes and remove supernatant.
- 7) Suspend cell pellet in 20mL of PBS and count PBMCs with a hemocytometer.
- 8) Determine what volume of the 20mL of cells to divide into a second conical tube for an unstained control. At this point one 50mL tube should have volume with approximately $20-35 \times 10^6$ cells for full staining and sorting. The second tube should have a smaller volume with approximately $5-10 \times 10^6$ cells to be used as an unstained control.
- 9) Spin at 1500 rpm for 5 mins and remove supernatant.
- 10) Resuspend each tube of cells in 200uL of cold FACS Wash and transfer into FACS tubes for staining.
 - a. Unstained tube – 5-10 million cells
 - b. Full stain – 20-35 million cells

ISOLATED PBMC HUMAN MONOCYTES

**All reagents, centrifuges and cells must be kept at 4C*

- 1) Add thawed cells to 15mL of cold RPMI Media.
- 2) Spin down at 1500rpm for 5 minutes and remove supernatant.
- 3) Suspend cell pellet in 2.5 mL of PBS and count PBMCs with a hemocytometer.
- 4) Spin down at 1500rpm for 5 minutes and remove supernatant.

- 5) Wash by adding PBS up to 2.5mL, count cells, then spin down and remove supernatant.
- 6) Resuspend each tube of cells in 200uL of cold FACS Wash and transfer into FACS tubes for staining. (Leave out a small volume (~20 μ L) in a separate FACS Tube as an unstained control)

STAINING MONOCYTES

- 1) Add 100uL of 3X surface stain to cells in FACS tube and incubate in the dark at 4°C for 30 minutes. (Stain is optimized for 10-35 million cells (ideally 20 million)/100uL of stain).

| Antigen | Fluorophore | Antibody Concentration in 3x stain | One Sample Full Stain |
|-----------------------|-------------|------------------------------------|-------------------------------|
| CD45 | BV510 | 1/50 | 2 μ L |
| CD16 | PECy7 | 1/50 | 2 μ L |
| CD14 | BV421 | 1/50 | 2 μ L |
| CD11b | APC | 1/25 | 4 μ L |
| HLA-DR | PerCpCy5.5 | 1/25 | 4 μ L |
| CD15 | AF700 | 1/50 | 2 μ L |
| CD3 | AF700 | 1/25 | 4 μ L |
| CD56 | AF700 | 1/25 | 4 μ L |
| CD19 | AF700 | 1/25 | 4 μ L |
| Total Antibody | | | 28.0 μL |
| FACS Wash | | | 72.0 μL |

- 2) Cool centrifuge to 4°C.
- 3) To wash, top up tube with 4mL with FACS Wash, spin down and remove supernatant.
- 4) Top up tube with 4mL of RBC Lysis Buffer and incubate in the dark at 4°C for 5 minutes.
- 5) Spin at 1500 rpm for 5 mins and remove supernatant.
- 6) Add 2mL of FACS Wash and count cells.
- 7) Spin at 1500 rpm for 5 mins and remove supernatant.

PREPARTION FOR CELL SORT

- 1) Prepare for cell sorting by resuspending in 1.0 mL of FACS Wash.
- 2) Prepare one tube with 20uL of fully stained cells as an 7AAD FMO control (add 220uL of FACS Wash)
- 3) Add 7AA-D at 1/250 (V/V) (4 μ L) to each sample. Mix with pipette.

- 4) Prepare one tube for cytometer set up (20ul of fully stained cells (with 7AAD) and 220 FACS Wash)
- 5) Prepare compensation controls. Note: 7AAD compensation is made with unstained cells rather than beads.
- 6) Set up voltages and compensation on sorter.
- 7) Filter cells through 0.45um mesh just prior to sorting.
- 8) Coat 2mL polypropylene FACS Tubes with FACS Wash (Cells will bust if hitting dry tubes).
- 9) Empty tubes before the cell sort.
- 10) Sort ~50 000 cells directly into 700uL QIAzol. Note that every 2000 cells are sorted in ~6uL of sorting buffer.
- 11) Verify viability and purity of cells post-sorting.
 - **Plan 1 hour of set up time on the cytometer
 - **Plan 1.5-2hours/sample of sort time.

ANTIBODY LIST

| Antibody | Fluorophore | Company | Clone | Catalog # |
|----------|-------------|---------------|----------|------------|
| CD45 | BV510 | BioLegend | HI30 | 304036 |
| CD16 | PE-Cy7 | eBioscience | eBioCB16 | 25-0168-42 |
| CD14 | BV421 | BioLegend | M5E2 | 301830 |
| CD11b | APC | BD Pharmingen | ICRF44 | 550019 |
| HLA-DR | PerCPCy5.5 | eBioScience | LN3 | 45-9956-42 |
| CD15 | AF700 | BioLegend | W6D3 | 323026 |
| CD3 | AF700 | BD Pharmingen | UCHT1 | 557943 |
| CD19 | AF700 | invitrogen | HIB19 | 56-0199-42 |
| CD56 | AF700 | BioLegend | HCD56 | 318316 |

SORTING STRATEGY

