



NF- κ B LUCIFERASE ASSAY

Sarah Chauvin: May 15th 2011

Updated by Charles Yin: July 11th, 2012

Bowdish Lab, McMaster University

Hamilton, ON, Canada

www.bowdish.ca

BACKGROUND

- The **luciferase reporter assay** is a technique in molecular biology with a diverse range of applications, including the characterization of signaling pathways. In the Bowdish Lab, we use the luciferase assay to measure NF- κ B levels in *in vitro* models of infection and inflammation.
- Luciferase is a family of enzymes produced naturally by a number of species, most notably those in the genus Lampyridae (firefly). These enzymes catalyze the oxidation of small, heterocyclic substrates called **luciferins**, resulting in the conversion of luciferin to an excited electronic state. When luciferin reverts back to ground state, luminescence results. This luminescence can then be detected by a luminometer.
- Molecular biologists take advantage of the luminescent property of luciferase by placing its gene under the control of promoters for a gene of interest. Subsequently, the transcriptional activation of this gene of interest can be monitored by luminescence.

NOTES

- This protocol is written for a procedure wherein HEK 293T cells are transiently transfected and then infected with *Streptococcus pneumoniae*. As such, it will be primarily of interest to Bowdish Lab members. **External users should modify this protocol where necessary.**
- **Plan for 5 consecutive days of lab work.** Day 1 to seed plates, Day 2 to transfect, Day 3 to infect and Day 5 to collect lysates and assay for luciferase activity.

EQUIPMENT

- Equipment:
 - o Plate reader **with an injection system** capable of reading luminescence
 - o **Opaque**, flat-bottom 96-well plates
 - o Multi-channel pipette capable of operating in the 10-50 μ l range (optional, but would make your life much easier)
- Materials:
 - o PEI transfection reagent
 - o Luminescent β -galactosidase Detection Kit II (Clontech Laboratories, Cat# 631712)
 - o Luciferase Assay Kit (Agilent Technologies, Cat# 219020)
 - o DMEM (all of: complete, serum-free, 1% FBS)

PROTOCOL

Part I: Seeding Cells

- Preparatory Work:
 - o Warm complete DMEM media and PBS to 37°C in water bath.
 - o Have ready a source of growing HEK 293T cells from which to seed.
 - o Prepare a BSL2 Biosafety Cabinet for work.
- 1. Remove media from growing HEK 293T cells.
- 2. Wash 1× with PBS.
- 3. Mechanically suspend cells in 20ml fresh media by washing bottom of flask repeatedly.
- 4. Assay cell count and viability. Typically, this is done through a Trypan blue stain.
- 5. Seed 6-well plates at a cell density of 3×10^5 cells per well in 3ml media.
 - a. **Note: Do not plate at a density lower than 3×10^5 cells, or they will not grow!**
- 6. Incubate overnight at 37°C. Cells should be 50-80% confluent after 24h.
- Concluding Work:
 - o Clean up BSL2 Biosafety Cabinet. Bleach and dispose of liquid waste.

Part II: Transfection

- Preparatory Work:
 - o Warm **serum-free** DMEM media and PBS to 37°C in water bath.
 - o Calculate the volumes of DNA needed for transfection. PEI transfection is optimal with a total of 2µg DNA. **Use extra to account for pipetting errors!**
 - o Prepare a BSL2 Biosafety Cabinet for work.
- 1. Add DNA to incomplete DMEM, using 100µl media per well to transfect (**and extra to account for pipetting error**).
- 2. Add 12µl 1ng/ml PEI per 100µl media. Vortex 15s. Incubate 10min at room temperature.
- 3. Add a further 600µl media per 100µl media.
- 4. Add 700µl media with DNA and PEI **evenly** and **drop-wise** per well.
- 5. Incubate 3h at 37°C.
- 6. Add 2.7ml complete DMEM media.
- 7. Incubate 24h at 37°C.
- Concluding Work:
 - o Clean up BSL2 Biosafety Cabinet. Bleach and dispose of liquid waste.

Part III: Infection

- Preparatory Work:
 - o Ensure that you have a heat-killed, lysozyme digested *S. pneumoniae* strain 1121 stock ready. **The concentration that our stocks are made at is 5×10^8 cells/ml.**
 - o Warm **1% FBS** DMEM media to 37°C in water bath.
 - o Calculate the volume of *S. pneumoniae* stock to use. **The optimal MOI is 25.**
 - o Prepare a BSL2 Biosafety Cabinet for work.
- 1. Add appropriate volume of *S. pneumoniae* stock to 1% FBS DMEM (3ml media per well).
- 2. Remove media from transfected HEK 293T cells.
- 3. Wash 1× with PBS, 1ml per well.

4. Add 3ml of media with *S. pneumoniae*.
5. Incubate 48h at 37°C.
- Concluding Work:
 - o Clean up BSL2 Biosafety Cabinet. Bleach and dispose of liquid waste.

Part IV: Collect Cell Lysates

- Preparatory Work:
 - o Prepare 1× lysis buffer from MilliQ H₂O 5× lysis buffer supplied by the Luciferase Assay kit.
 - o Prepare a bucket of ice. Bring with you to BSL2 Safety Cabinet.
 - o Prepare BSL2 Biosafety Cabinet for work.
- 1. Remove media from HEK 293T cells.
- 2. Wash 2× with PBS, 1ml per well.
- 3. Cover cells with 300µl 1× lysis buffer.
- 4. Incubate 15min at room temperature. Swirl every ~3min.
- 5. Scrape cells into suspension in lysis buffer. Transfer to a 1.5ml Eppendorf tube on ice.
- 6. Vortex tubes 15s.
- 7. Centrifuge 2min at 12,000×g and 4°C.
- 8. Transfer supernatant to fresh 1.5ml Eppendorf tube on ice. Dispose of pellet.
- 9. Assay immediately (**Part V**) or store at -80°C.
- Concluding Work:
 - o Clean BSL2 Biosafety Cabinet. Bleach and dispose of liquid waste.

Part V: Measuring Luminescence

- Preparatory Work:
 - o Prepare luciferase substrate-assay by adding the entire assay buffer into the vial containing the lyophilized luciferase substrate (provided by Luciferase Assay kit). Mix well and divide into 1ml aliquots. Store at -80°C. **You should make at least 1.3ml more than necessary to account for dead volume in the plate reader injection system.**
 - o Prepare β-galactosidase reaction substrate by warming buffer and substrate to room temperature and adding 1µl of substrate per 49µl buffer.
 - o **Book the plate reader with injection system well beforehand!**
- 1. Using an opaque 96-well plate, transfer 20µl cell lysate per well in triplicate to test for luciferase activity.
- 2. Transfer 10µl cell lysate per well in triplicate to test for β-galactosidase activity. **This is separate from the lysates for luciferase activity.**
- 3. Add 50µl of β-galactosidase reaction mixture to each well used to test for β-galactosidase activity.
- 4. Incubate plate 30min at room temperature covered with aluminium foil.
- 5. Meanwhile, transfer enough luciferase reaction mixture (you will be using 50µl per well) to a 15ml Falcon tube (**adding enough to account for injector dead volume**). Cover with aluminium foil.
- 6. Measure luminescence of wells on a plate reader. For luciferase measurements, set injector to inject 50µl reaction mixture **immediately** prior to luminescence measurement.
- Concluding Work:
 - o **Note: β-galactosidase is used as a transfection control and luciferase luminescence measurements are normalized to β-galactosidase luminescence measurements.**